

Product Datasheet

CD107b antibody (FITC) (orb433886)

Description	Mouse monoclonal antibody against CD107b conjugated to FITC
Species/Host	Mouse
Reactivity	Canine
Conjugation	FITC
Tested Applications	FC
Immunogen	MDCK (Madin-Darby Canine Kidney) cells
Target	CD107b
Preservatives	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.1mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.