

Product Datasheet

CD134 antibody (FITC) (orb433893)

Description	Mouse monoclonal antibody against CD134 conjugated to FITC
Species/Host	Mouse
Reactivity	Feline
Conjugation	FITC
Tested Applications	FC
Immunogen	CHO-derived feline CD134-Fc fusion protein.
Target	CD134
Preservatives	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.1mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.