

Product Datasheet

HERPES SIMPLEX VIRUS 1/2 gD antibody (orb433900)

Description	Mouse monoclonal antibody to HERPES SIMPLEX VIRUS 1/2 gD
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF
Target	HERPES SIMPLEX VIRUS 1/2 gD
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	0.1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.