

Product Datasheet

Coronavirus antibody (orb433926)

Description	Mouse monoclonal antibody against CORONAVIRUS conjugated to Biotin
Species/Host	Mouse
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA
Target	CORONAVIRUS
Preservatives	Purified IgG conjugated to Biotin - liquid
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.