

Product Datasheet

FK-506 antibody (orb433954)

Description	Mouse monoclonal antibody to FK-506
Species/Host	Mouse
Reactivity	Chemical
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	FK-506-BSA conjugate.
Target	FK-506
Preservatives	Purified IgM - liquid
Form/Appearance	Purified IgM - liquid
Concentration	1.0mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgM
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.