

Product Datasheet

RESPIRATORY SYNCYTIAL VIRUS FUSION PROTEIN antibody (orb433998)

Description	Mouse monoclonal antibody to RESPIRATORY SYNCYTIAL VIRUS FUSION PROTEIN
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF
Immunogen	Fusion protein of the human Respiratory Syncytial Virus.
Target	RESPIRATORY SYNCYTIAL VIRUS FUSION PROTEIN
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	0.1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.