

Product Datasheet

INFLUENZA A NUCLEOPROTEIN antibody (orb434133)

Description	Mouse monoclonal antibody to INFLUENZA A NUCLEOPROTEIN
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IH-Fr
Immunogen	Purified influenza virus type A (H1N1).
Target	INFLUENZA A NUCLEOPROTEIN
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.