

## Product Datasheet

### HEPATITIS C CORE ANTIGEN antibody (orb434173)

<b>Description</b>	Mouse monoclonal antibody to HEPATITIS C CORE ANTIGEN
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, ICC, IF, WB
<b>Immunogen</b>	Recombinant full-length core protein.
<b>Target</b>	HEPATITIS C CORE ANTIGEN
<b>Preservatives</b>	Purified IgG - liquid
<b>Form/Appearance</b>	Purified IgG - liquid
<b>Concentration</b>	0.1mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Purified
<b>Expiration Date</b>	12 months from date of receipt.