

## Product Datasheet

# PAPILLOMAVIRUS 16 ONCOPROTEIN E7 antibody (orb434198)

<b>Description</b>	Mouse monoclonal antibody to PAPILLOMAVIRUS 16 ONCOPROTEIN E7
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	Oncoprotein E7 from human papillomavirus 16.
<b>Target</b>	PAPILLOMAVIRUS 16 ONCOPROTEIN E7
<b>Preservatives</b>	Purified IgG - liquid
<b>Form/Appearance</b>	Purified IgG - liquid
<b>Concentration</b>	1.0mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Purified
<b>Uniprot ID</b>	<b>P03129</b>
<b>Expiration Date</b>	12 months from date of receipt.