

Product Datasheet

PAPILLOMAVIRUS 16 ONCOPROTEIN E7 antibody (orb434199)

Description	Mouse monoclonal antibody to PAPILLOMAVIRUS 16 ONCOPROTEIN E7
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Oncoprotein E7 from human papillomavirus 16,
Target	PAPILLOMAVIRUS 16 ONCOPROTEIN E7
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P03129
Expiration Date	12 months from date of receipt.