

Product Datasheet

S-TAG antibody (orb434425)

Description	Mouse monoclonal antibody to S-TAG
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FACS, Purification
Immunogen	S-tag recombinant protein.
Target	S-TAG
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	0.5mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.