

Product Datasheet

METAPNEUMOVIRUS antibody (orb434429)

| | |
|----------------------------|---|
| Description | Mouse monoclonal antibody to METAPNEUMOVIRUS |
| Species/Host | Mouse |
| Reactivity | Virus |
| Conjugation | Unconjugated |
| Tested Applications | ICC, IF |
| Immunogen | Sub-type A hMPV virus isolate NCL03-4/145. |
| Target | METAPNEUMOVIRUS |
| Preservatives | Purified IgG - liquid |
| Form/Appearance | Purified IgG - liquid |
| Concentration | 1.0mg/ml |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG2a |
| Clonality | Monoclonal |
| Purity | Purified |
| Expiration Date | 12 months from date of receipt. |