

Product Datasheet

DYKDDDDK TAG antibody (orb434542)

Description	Rat monoclonal antibody to DYKDDDDK TAG
Species/Host	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FACS, ICC, IF, IHC-P, IP, WB
Target	DYKDDDDK TAG
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.