

## Product Datasheet

# RESPIRATORY SYNCYTIAL VIRUS FUSION PROTEIN antibody (orb434582)

<b>Description</b>	Mouse monoclonal antibody to RESPIRATORY SYNCYTIAL VIRUS FUSION PROTEIN
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, FACS, ICC, IF, IP, WB
<b>Immunogen</b>	Bovine RSV strains: 127, SNK and 9007. Human RSV strains: Long, Randall, 8/60, and A/2.
<b>Target</b>	RESPIRATORY SYNCYTIAL VIRUS FUSION PROTEIN
<b>Preservatives</b>	Purified IgG - liquid
<b>Form/Appearance</b>	Purified IgG - liquid
<b>Concentration</b>	1.0 mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Purified
<b>Expiration Date</b>	12 months from date of receipt.