

Product Datasheet

RESPIRATORY SYNCYTIAL VIRUS NUCLEOPROTEIN antibody (orb434583)

Description	Mouse monoclonal antibody to RESPIRATORY SYNCYTIAL VIRUS NUCLEOPROTEIN
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr
Immunogen	Bovine RSV strains: 127, SNK and 9007. Human RSV strains: A/2, Long, Randall and 8/60.
Target	RESPIRATORY SYNCYTIAL VIRUS NUCLEOPROTEIN
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.