

Product Datasheet

MACROPHAGES/MONOCYTES antibody (FITC) (orb434632)

Description	Rat monoclonal antibody against MACROPHAGES/MONOCYTES conjugated to FITC
Species/Host	Rat
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	Mouse lymph node stroma.
Target	MACROPHAGES/MONOCYTES
Preservatives	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.