

Product Datasheet

HISTIDINE TAG antibody (HRP) (orb434953)

Description	Mouse monoclonal antibody against HISTIDINE TAG conjugated to HRP
Species/Host	Mouse
Conjugation	HRP
Tested Applications	ELISA, WB
Target	HISTIDINE TAG
Preservatives	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Form/Appearance	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.