

## Product Datasheet

### CD134 antibody (FITC) (orb435082)

<b>Description</b>	Mouse monoclonal antibody against CD134 conjugated to FITC
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Rat
<b>Conjugation</b>	FITC
<b>Tested Applications</b>	FC
<b>Immunogen</b>	Phytohemagglutinin (PHA) activated Rat Lymph Node Cells.
<b>Target</b>	CD134
<b>Preservatives</b>	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
<b>Form/Appearance</b>	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
<b>Concentration</b>	0.1 mg/ml
<b>Storage</b>	+4°C, -20°C if preferred
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<b>P15725</b>
<b>Expiration Date</b>	12 months from date of receipt.