

Product Datasheet

MACROPHAGES/MONOCYTES antibody (FITC) (orb435213)

| | |
|----------------------------|--|
| Description | Mouse monoclonal antibody against MACROPHAGES/MONOCYTES conjugated to FITC |
| Species/Host | Mouse |
| Reactivity | Human |
| Conjugation | FITC |
| Tested Applications | FC |
| Immunogen | Human monocytes. |
| Target | MACROPHAGES/MONOCYTES/GRANULOCYTES |
| Preservatives | Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid |
| Form/Appearance | Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid |
| Concentration | 0.1 mg/ml |
| Storage | +4°C, -20°C if preferred |
| Note | For research use only |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Uniprot ID | P06702 |
| Expiration Date | 12 months from date of receipt. |