

Product Datasheet

MACROPHAGES/MONOCYTES antibody (FITC) (orb435214)

Description	Mouse monoclonal antibody against MACROPHAGES/MONOCYTES conjugated to FITC
Species/Host	Mouse
Reactivity	Human
Conjugation	FITC
Tested Applications	FC
Immunogen	Human monocytes.
Target	MACROPHAGES/MONOCYTES/GRANULOCYTES
Preservatives	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P06702
Expiration Date	12 months from date of receipt.