

Product Datasheet

CYTOMEGALOVIRUS IMMEDIATE/EARLY ANTIGEN antibody (orb435332)

Description	Mouse monoclonal antibody to CYTOMEGALOVIRUS IMMEDIATE/EARLY ANTIGEN
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF
Target	CYTOMEGALOVIRUS IMMEDIATE/EARLY ANTIGEN
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	0.1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.