

Product Datasheet

PARAINFLUENZA VIRUS 1 antibody (orb435362)

Description	Goat polyclonal antibody to PARAINFLUENZA VIRUS 1
Species/Host	Goat
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Parainfluenza type 1 - Cantell strain.
Target	PARAINFLUENZA VIRUS 1
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.