

Product Datasheet

INFLUENZA A antibody (orb435398)

Description	Mouse monoclonal antibody to INFLUENZA A
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Influenza A blend.
Target	INFLUENZA A
Preservatives	Ascites - liquid
Form/Appearance	Ascites - liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Cocktail
Clonality	Monoclonal
Purity	Ascites
Expiration Date	12 months from date of receipt.