

Product Datasheet

INFLUENZA A H1 ANTIGEN antibody (orb435404)

Description	Mouse monoclonal antibody to INFLUENZA A H1 ANTIGEN
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Influenza blend.
Target	INFLUENZA A H1 ANTIGEN
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.