

Product Datasheet

INFLUENZA HEMAGGLUTININ antibody (orb435413)

Description	Chicken polyclonal antibody to INFLUENZA HEMAGGLUTININ
Species/Host	Gallus
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	IH-Fr, IP, WB
Immunogen	Synthetic peptide YPYDVPDYA from influenza hemagglutinin (HA).
Target	INFLUENZA HEMAGGLUTININ
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	0.1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgY
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.