

Product Datasheet

ADENOVIRUS HEXON antibody (FITC) (orb435436)

Description	Mouse monoclonal antibody against ADENOVIRUS HEXON conjugated to FITC
Species/Host	Mouse
Reactivity	Virus
Conjugation	FITC
Tested Applications	ELISA, ICC, IF
Target	ADENOVIRUS HEXON
Preservatives	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.