

## Product Datasheet

### Viral ADENOVIRUS protein (orb435586)

<b>Description</b>	Native viral ADENOVIRUS protein
<b>Reactivity</b>	Virus
<b>Tested Applications</b>	ELISA
<b>Target</b>	ADENOVIRUS
<b>Preservatives</b>	Semi pure suspension of Adenoid 6 strain, grade 2 - liquid
<b>Form/Appearance</b>	Semi pure suspension of Adenoid 6 strain, grade 2 - liquid
<b>Concentration</b>	3.24 mg/ml
<b>Storage</b>	-70°C (ship frozen)
<b>Note</b>	For research use only
<b>Purity</b>	Purified
<b>Expiration Date</b>	6 months from date of receipt.