

## Product Datasheet

### Human C-REACTIVE PROTEIN protein (orb435594)

<b>Description</b>	Native human C-REACTIVE PROTEIN protein
<b>Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Target</b>	C-REACTIVE PROTEIN
<b>Preservatives</b>	Purified protein from human plasma - liquid
<b>Form/Appearance</b>	Purified protein from human plasma - liquid
<b>Concentration</b>	1.0 mg/ml
<b>Storage</b>	+4°C, -20°C if preferred
<b>Note</b>	For research use only
<b>Purity</b>	Purified
<b>Uniprot ID</b>	<b>P02741</b>
<b>Expiration Date</b>	6 months from date of receipt.