

## Product Datasheet

### Human FERRITIN protein (orb435614)

<b>Description</b>	Native human FERRITIN protein
<b>Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Target</b>	FERRITIN
<b>Preservatives</b>	Purified protein from human liver - liquid
<b>Form/Appearance</b>	Purified protein from human liver - liquid
<b>Concentration</b>	3.0 mg/ml
<b>Storage</b>	+4°C DO NOT FREEZE
<b>Note</b>	For research use only
<b>Purity</b>	Purified
<b>Uniprot ID</b>	<b>P02792, P02794</b>
<b>Expiration Date</b>	6 months from date of receipt.