

Product Datasheet

Human D-DIMER protein (orb435616)

Description	Native human D-DIMER protein
Reactivity	Human
Tested Applications	ELISA
Target	D-DIMER
Preservatives	Purified protein from human plasma - liquid
Form/Appearance	Purified protein from human plasma - liquid
Concentration	1.0mg/ml
Storage	-20°C only (ship frozen)
Note	For research use only
Purity	Purified
Uniprot ID	P02671, P02679, P02675
Expiration Date	6 months from date of receipt.