

Product Datasheet

Human FIBRIN DEGRADATION PRODUCT E protein (orb435618)

Description	Native human FIBRIN DEGRADATION PRODUCT E protein
Reactivity	Human
Tested Applications	ELISA
Target	FIBRIN DEGRADATION PRODUCT E
Preservatives	Purified protein from human serum - lyophilized
Form/Appearance	Purified protein from human serum - lyophilized
Concentration	0.1 mg/ml
Storage	-20°C after reconstitution
Note	For research use only
Purity	Purified
Uniprot ID	P02671, P02679, P02675
Expiration Date	6 months from date of receipt.