

Product Datasheet

Bovine HYALURONIDASE protein (orb435625)

Description	Native bovine HYALURONIDASE protein
Reactivity	Bovine
Tested Applications	ELISA
Target	HYALURONIDASE
Preservatives	Purified enzyme from bovine testes - liquid
Form/Appearance	Purified enzyme from bovine testes - liquid
Concentration	0.2 g/ml
Storage	-20°C only (ship frozen)
Note	For research use only
Purity	Purified
Uniprot ID	Q8SQG8 , Q5E985
Expiration Date	6 months from date of receipt.