

Product Datasheet

Human ALPHA 1 MICROGLOBULIN protein (orb435637)

Description	Native human ALPHA 1 MICROGLOBULIN protein
Reactivity	Human
Tested Applications	ELISA
Target	ALPHA 1 MICROGLOBULIN
Preservatives	Purified protein from the urine of patients with chronic renal tubular proteinuria - lyophilized
Form/Appearance	Purified protein from the urine of patients with chronic renal tubular proteinuria - lyophilized
Concentration	1.0mg/ml
Storage	-20°C after reconstitution
Note	For research use only
Purity	Purified
Uniprot ID	P02760
Expiration Date	6 months from date of receipt.