

Product Datasheet

Rabbit PHOSPHORYLASE B protein (orb435645)

Description	Native rabbit PHOSPHORYLASE B protein
Reactivity	Rabbit
Tested Applications	ELISA
Target	PHOSPHORYLASE B
Preservatives	Purified enzyme from rabbit skeletal muscle - liquid
Form/Appearance	Purified enzyme from rabbit skeletal muscle - liquid
Concentration	5.0 mg/ml
Storage	-20°C only (ship frozen)
Note	For research use only
Purity	Purified
Expiration Date	6 months from date of receipt.