

Product Datasheet

Human PROSTATE SPECIFIC ANTIGEN protein (orb435649)

Description	Native human PROSTATE SPECIFIC ANTIGEN protein
Reactivity	Human
Tested Applications	ELISA, WB
Target	PROSTATE SPECIFIC ANTIGEN
Preservatives	Purified protein from human seminal fluid - liquid
Form/Appearance	Purified protein from human seminal fluid - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Purity	Purified
Uniprot ID	P07288
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)