

Product Datasheet

Human PROSTATE SPECIFIC ANTIGEN protein (orb435650)

| | |
|----------------------------|--|
| Description | Native human PROSTATE SPECIFIC ANTIGEN protein |
| Reactivity | Human |
| Tested Applications | ELISA |
| Target | PROSTATE SPECIFIC ANTIGEN |
| Preservatives | Purified protein from human seminal fluid - liquid |
| Form/Appearance | Purified protein from human seminal fluid - liquid |
| Concentration | 1.0 mg/ml |
| Storage | +4°C, -20°C if preferred |
| Note | For research use only |
| Purity | Purified |
| Uniprot ID | P07288 |
| Expiration Date | 6 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)