

Product Datasheet

Human RETINOL BINDING PROTEIN protein (orb435656)

Description	Native human RETINOL BINDING PROTEIN protein
Reactivity	Human
Tested Applications	ELISA
Target	RETINOL BINDING PROTEIN
Preservatives	Purified protein from the urine of patients with chronic renal tubular proteinuria - lyophilized
Form/Appearance	Purified protein from the urine of patients with chronic renal tubular proteinuria - lyophilized
Concentration	1.0 mg/ml
Storage	-20°C after reconstitution
Note	For research use only
Purity	Purified
Uniprot ID	P02753
Expiration Date	6 months from date of receipt.