

Product Datasheet

Human TENASCIN protein (orb435659)

Description	Native human TENASCIN protein
Reactivity	Human
Tested Applications	ELISA
Target	TENASCIN
Preservatives	Purified protein from tissue culture supernatant - liquid
Form/Appearance	Purified protein from tissue culture supernatant - liquid
Concentration	0.1mg/ml
Storage	-70°C (ship frozen)
Note	For research use only
Purity	Purified
Uniprot ID	P24821
Expiration Date	6 months from date of receipt.