

Product Datasheet

Human THYROGLOBULIN protein (orb435660)

Description	Native human THYROGLOBULIN protein
Reactivity	Human
Tested Applications	ELISA
Target	THYROGLOBULIN
Preservatives	Purified protein from human thyroid glands - liquid
Form/Appearance	Purified protein from human thyroid glands - liquid
Concentration	1.0 mg/ml
Storage	-20°C only (ship frozen)
Note	For research use only
Purity	Purified
Uniprot ID	P01266
Expiration Date	6 months from date of receipt.