

Product Datasheet

Chicken IL6 protein (orb435709)

Description	Recombinant of chicken IL6 protein
Reactivity	Gallus
Tested Applications	FA
Target	IL-6
Preservatives	Purified His-tagged recombinant protein expressed in E.coli - liquid
Form/Appearance	Purified His-tagged recombinant protein expressed in E.coli - liquid
Concentration	500 µg/ml
Storage	-20°C only (ship frozen)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Uniprot ID	Q90Y10
Expiration Date	6 months from date of receipt.