

Product Datasheet

Bovine IFN GAMMA protein (orb435717)

Description	Recombinant of bovine IFN GAMMA protein
Reactivity	Bovine
Tested Applications	ELISA, FA
Target	IFN GAMMA
Preservatives	Purified recombinant protein - lyophilized
Form/Appearance	Purified recombinant protein - lyophilized
Storage	-20°C after reconstitution
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Uniprot ID	P07353
Expiration Date	6 months from date of receipt.