

Product Datasheet

Bovine IFN ALPHA A protein (orb435726)

Description	Recombinant of bovine IFN ALPHA A protein
Reactivity	Bovine
Tested Applications	ELISA, FA, WB
Target	IFN ALPHA A
Preservatives	Purified recombinant protein - lyophilized
Form/Appearance	Purified recombinant protein - lyophilized
Storage	-20°C only (ship +4°C)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Uniprot ID	P05007
Expiration Date	6 months from date of receipt.