

Product Datasheet

Bovine IL6 protein (orb435728)

Description	Recombinant of bovine IL6 protein
Reactivity	Bovine
Tested Applications	ELISA, FA, WB
Target	IL-6
Preservatives	Purified recombinant protein - lyophilized
Form/Appearance	Purified recombinant protein - lyophilized
Storage	-20°C only (ship +4°C)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Uniprot ID	P26892
Expiration Date	6 months from date of receipt.