

Product Datasheet

Bovine CXCL9 protein (orb435729)

Description	Recombinant of bovine CXCL9 protein
Reactivity	Bovine
Tested Applications	ELISA
Target	CXCL9
Preservatives	Purified recombinant protein - lyophilized
Form/Appearance	Purified recombinant protein - lyophilized
Storage	-20°C only (ship +4°C)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Uniprot ID	A9QWP9
Expiration Date	6 months from date of receipt.