

Product Datasheet

Bovine IL1 ALPHA protein (orb435736)

Description	Recombinant of bovine IL1 ALPHA protein
Reactivity	Bovine
Tested Applications	ELISA, FA
Target	IL-1 ALPHA
Preservatives	Purified recombinant protein - lyophilized
Form/Appearance	Purified recombinant protein - lyophilized
Storage	-20°C only (ship +4°C)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Uniprot ID	P08831
Expiration Date	6 months from date of receipt.