

## Product Datasheet

### Bovine IFN BETA protein (orb435737)

<b>Description</b>	Recombinant of bovine IFN BETA protein
<b>Reactivity</b>	Bovine
<b>Tested Applications</b>	ELISA, WB
<b>Target</b>	IFN BETA
<b>Preservatives</b>	Purified recombinant protein - lyophilized
<b>Form/Appearance</b>	Purified recombinant protein - lyophilized
<b>Storage</b>	-20°C only (ship +4°C)
<b>Note</b>	For research use only
<b>Clonality</b>	Recombinant
<b>Purity</b>	Rec. Protein
<b>Expiration Date</b>	6 months from date of receipt.