

Product Datasheet

Human KAPPA LIGHT CHAIN protein (orb435888)

Description	Native human KAPPA LIGHT CHAIN protein
Reactivity	Human
Tested Applications	ELISA
Target	KAPPA LIGHT CHAIN
Preservatives	Purified human kappa light chain - liquid
Form/Appearance	Purified human kappa light chain - liquid
Concentration	5.0 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Purity	Preservative Free
Uniprot ID	P01834
Expiration Date	6 months from date of receipt.