

## Product Datasheet

### DYKDDDDK TAG antibody (HRP) (orb436316)

<b>Description</b>	Rabbit polyclonal antibody against DYKDDDDK TAG conjugated to HRP
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	Keyhole limpet haemocyanin conjugated Enterokinase cleavage site peptide DYKDDDDK. Residues of cysteine and glycine were used to facilitate coupling at the C-terminal end.
<b>Target</b>	DYKDDDDK TAG
<b>Preservatives</b>	Purified IgG conjugated to Horseradish Peroxidase (HRP) - lyophilized
<b>Form/Appearance</b>	Purified IgG conjugated to Horseradish Peroxidase (HRP) - lyophilized
<b>Concentration</b>	1mg/ml
<b>Storage</b>	-20°C after reconstitution
<b>Note</b>	For research use only
<b>Isotype</b>	Polyclonal IgG
<b>Clonality</b>	Polyclonal
<b>Expiration Date</b>	12 months from date of receipt.