

## Product Datasheet

### HISTIDINE TAG antibody (FITC) (orb436319)

<b>Description</b>	Goat polyclonal antibody against HISTIDINE TAG conjugated to FITC
<b>Species/Host</b>	Goat
<b>Conjugation</b>	FITC
<b>Tested Applications</b>	ICC, IF
<b>Immunogen</b>	Keyhole limpet haemocyanin (KLH) conjugated H-H-H-H-H
<b>Target</b>	HISTIDINE TAG
<b>Preservatives</b>	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
<b>Form/Appearance</b>	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
<b>Concentration</b>	1.0mg/ml
<b>Storage</b>	+4°C DO NOT FREEZE
<b>Note</b>	For research use only
<b>Isotype</b>	Polyclonal IgG
<b>Clonality</b>	Polyclonal
<b>Expiration Date</b>	12 months from date of receipt.