

Product Datasheet

V5-TAG antibody (HRP) (orb436354)

Description	Mouse monoclonal antibody against V5-TAG conjugated to HRP
Species/Host	Mouse
Reactivity	Virus
Conjugation	HRP
Tested Applications	ELISA, WB
Immunogen	Paramyxovirus Simian-Virus 5 (SV5)
Target	V5-TAG
Preservatives	Purified IgG conjugated to horseradish peroxidase (HRP) - liquid
Form/Appearance	Purified IgG conjugated to horseradish peroxidase (HRP) - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P11207
Expiration Date	12 months from date of receipt.